



Release Notes

NX Tester - Release 13 Serial Port Upgrade

October, 2009

Component Versions

- NX Tester Firmware Version: 2.13.149

WARNING:

- Do not run this serial port upgrade on any tester that has a version newer than Release 12 !!!

If Current Version is Release 12 (April 2009) or older, then the NX Tester must be upgraded to Release 13 via the serial port:

- 1) Connect your NX Tester Serial Port (Serial 1) to your PC using the Dynalab download cable (included with the NX Editor)
- 2) Make sure the NX Tester is powered on and at the main menu
- 3) On your PC, open the zip file, and execute "NX_Tester_2_13_150.exe"
- 4) Follow the on-screen instructions
- 5) Then upgrade to the latest release via the memory card